

SAFETY NOTICE

SIT/STAND ROLLING WALKER ITEM# 66838

There is the possibility that a walker wheel on this product may come off following repeated use without tightening of the wheel bolt, presenting a hazard for the user.



Consumers should **STOP USING THESE WALKERS IMMEDIATELY**. Bring the walker along with a purchase receipt to your local Harbor Freight Tools' Retail Store for a refund. Call Harbor Freight Tools at **(800) 444-3353** or send an email to recall@harborfreight.com Monday through Friday between 8:00 am and 4:30 pm PDT to learn how to return the walker for a refund from Harbor Freight Tools.

Post until October 1, 2014

This Recall is being conducted following notice to the Food and Drug Administration